

# EXSUDEX®

wounddrainage

THE **MEDICAL &  
WOUNDCARE**  
COMPANY



## Indication for use

- Pressure ulcers
- Chronical wounds
- Ulcer: venous/arterial
- Diabetic ulser
- Infected wounds
- Acute & Traumatic wounds
- Partial thickness burns
- Dehisced incisions
- Meshed grafts
- Flaps

## Positively inexpensive

- Lower rentals
- Less frequent dressing changes
- Disposable canister
- Cost-effective consumables, either gauze or foam

## Positively simpler to use


- Simple setup
- Easy to use controls
- Fast and easy dressing changes
- Intermittend function
- Tube blocking alarm
- Standby alarm
- Menu bar with display of settings

## Positively portable


- Main or battery operation
- 12 hours duration

## Contra indications and precautions

- Malignancy in the wound.  
*However can be used for wound drainage in terminally ill patients.*
- Fistulae of unknown origin.  
*A non-enteric fistula or an unexplored fistula to organs or body cavities.*
- Necrotic tissue with eschar present.  
*Wound should be debrided to remove majority of necrotic tissue and eschar.*
- Untreated Osteomyelitis.  
*Necrotic bone should be removed  
Antibiotics should be administered.*
- Exposed blood vessels or organs.  
*Ensure all vessels or organs are protected with overlying fascia, tissue or other protective barrier.*
- Precautions should be taken with:  
*Patients on Anticoagulants.  
Patients with difficult wound homeostatis.  
Patients with active bleeding.*



## Technical specifications

Pressure range:	-50 mm Hg to 200 mmHg
Power Supply:	100-240 VAC 50-60 Hz, 350-150 mA 5V DC, 2.4A 
Max. use:	10 watt
Dimensions:	HxWxD(mm) 255x145x125 (In combination with canister)
Weight:	1.1 kg
Type (IEC 60601-1):	BF
Classification (IEC 60601-1):	II
Classification (93/42/EEC):	Ila



Produced in the Netherlands by MediTop BV

THE **MEDICAL & WOUND CARE** COMPANY

The Medical & Woundcare Company BV  
Vlasakker 22  
3417 XT Montfoort, The Netherlands  
Phone +31 348 47 58 49  
Fax +31 348 47 48 78  
E-mail: info@tmwc.eu  
www.themedicalwoundcarecompany.