

**EXSUDEX<sup>®</sup>**  
wounddrainage

THE **MEDICAL &  
WOUNDCARE**  
COMPANY

# EXSUDEX<sup>®</sup> XS



## Indication for use

- Pressure ulcers
- Chronical wounds
- Ulcer: venous/arterial
- Diabetic ulser
- Infected wounds
- Acute & Traumatic wounds
- Partial thickness burns
- Dehisced incisions
- Meshed grafts
- Flaps

## Positively inexpensive

- Lower rentals
- Less frequent dressing changes
- Disposable canister
- Cost-effective consumables, either gauze or foam

## Positively simpler to use

- Simple setup
- Easy to use controls
- Fast and easy dressing changes
- Intermittend function
- Tube blocking alarm
- Standby alarm
- Menu bar with display of settings

## Positively portable

- Main or battery operation
- 12 hours duration

## Contra indications and precautions

- Malignancy in the wound.

However can be used for wound drainage in terminally ill patients.

- Fistulae of unknown origin.

A non-enteric fistula or an unexplored fistula to organs or body cavities.

- Necrotic tissue with eschar present.

Wound should be debrided to remove majority of necrotic tissue and eschar.

- Untreated Osteomyelitis.

Necrotic bone should be removed

Antibiotics should be administered.

- Exposed blood vessels or organs.



Ensure all vessels or organs are protected with overlying fascia, tissue or other protective barrier.

- Precautions should be taken with:

Patients on Anticoagulants.

Patients with difficult wound homeostatis.

Patients with active bleeding.



## Technical specifications

Pressure range:	-10 mm Hg to 200 mmHg
Power Supply:	100-240 VAC 50-60 Hz
Internal Voltage	5V DC, 2A
Max. use:	10 watt
Dimensions:	HxWxD(mm) 148x151x94.5
Weight:	800 gram with empty canister 1125 gram with full canister

Type (IEC 60601-1): BF

Classification (IEC 60601-1): II

Classification (93/42/EEC): IIa

Produced in the Netherlands by MediTop BV

THE **MEDICAL &**  
**WOUND CARE**  
COMPANY

The Medical & Woundcare Company BV

Vlasakker 22

3417 XT Montfoort, The Netherlands

Phone +31 348 47 58 49

Fax +31 348 47 48 78

E-mail: info@tmwc.eu

www.themedicalwoundcarecompany.com

